



Applications of Two-Phase (Liquid/Gas) Xenon Gamma Cameras to the Detection of Special Nuclear Material and PET Medical Imaging

Nicholas E. Destefano
University of Connecticut
Yale University

- Noble Liquid Detector
- Applications of Xenon TPC Detectors
- The PIXeY Detector

DOE grant Number DE-SC0004599

ANS&T Washington, D.C. August 22, 2011

Applications of Two Phase (Liquid/Gas) Xenon Gamma Cameras to the Detection of Special Nuclear Material and PET Medical Imaging

Nicholas E. Destefano, Moshe Gai, Aaron G. Swindell, Alexander H. Young
University of Connecticut

Daniel N. McKinsey, Ethan P. Bernard, Sidney B. Cahn, Blair Edwards,
Christopher J. Kachulis, Nicole A. Larsen, Alexey Lyashenko, James A. Nikkel,
Christopher G. Wahl
Yale University



Detector Materials

	Energy resolution (σ @ 662 keV)	Position resolution (cm)	Density (g/cc)	Cost (\$/kg)	Detector element thickness (cm)	n/gamma discrimination
Plastic scintillator	Poor	~10	1	~30	~20	No
Liquid scintillator	17%	~10	0.9	~30	~20	Yes
Sodium Iodide	5%	~ few cm	3.7	~800	~10	No
CZT	0.33%	~0.1	5.8	~50,000	~1	No
Germanium	0.17%	~0.3	5.3	~20,000	~5	No
Liquid xenon	2.5-5%	~0.1	3	~1,000	~20	Yes


Liquid xenon combines:

- Energy resolution between CZT and sodium iodide
- Cost of sodium iodide
- Position resolution of CZT
- Scalability of organic scintillator


Additional Properties:

- 6-cm attenuation length at 1 MeV
- ~25 ns response
- 13.6 eV/e⁻

Xe Charge and Light Production




Two-phase Xe Detectors: How They Work




Project Goals

- Development and commissioning of cryogenics
- Monte Carlo simulations of Compton imaging
- Demonstration of superior energy resolution in dual phase operation




Noble Liquid Detector Applications


- Fundamental Physics
- Passive/Active Cargo Inspection of Special Nuclear Materials (SNM) With Compton Imaging
- Three Photon PET





Double Beta Decay



Two neutrino double beta decay ($2\nu\beta\beta$)





Neutrinoless double beta decay ($0\nu\beta\beta$)





Compton Imaging

$$E'_\gamma = \frac{E_\gamma \cdot m_e c^2}{E_\gamma \cdot (1 - \cos\varphi) + m_e c^2}$$

$$\varphi = \arccos \left[1 - m_e c^2 \cdot \left(\frac{1}{E_\gamma - E_1} - \frac{1}{E_\gamma} \right) \right]$$









Simulation of Three Photon PET

(With a Compton Gamma-Camera and a 3γ emitter)





- 3 photon PET: ~ 2.2 -mm spatial resolution
- Alternative to MicroPET with comparable resolution





voxel: $2 \times 2 \times 2$ mm³

The XENON10 Detector



Xenon10



Anti-correlation considerably improves energy resolution.

arXiv: 1001.2834

Nicholas E. Destefano, UConn and Yale University. ANS&T Washington D.C. August 22, 2011

PIXeY – Particle Identification in Xenon at Yale



- Separate drift and proportional scintillation fields
- Optimize energy resolution



Photomultiplier Tubes





- 14 two-inch Hamamatsu R8778 PMTs
- Quantum efficiency $\sim 35\%$




Mesh Grids




- Grids establish electric field, as well as protect PMTs
- Optical Transparency: 88%
- Completed construction, currently installed.



Parallel Single Wire Grids





- Higher strength for given optical transparency
 - Less grid deflection
 - Greater field uniformity
 - Better energy resolution
- Can be custom made for each application




Low Noise Wire Readout Charge Preamplifiers


- Read individual wires in LXe
- With 3 mm wire pitch, expect < 1 mm position resolution
- Gain: $5.33 \mu\text{V}/e^-$
- Noise: $120 e^-$
- Preamplifiers planned to be mounted inside detector, at LXe temperature
- Design used for GERDA germanium double beta decay experiment
- **Testing underway**




Monte Carlo Studies With GEANT4



- Ten detector geometries studied, to optimize efficiency and resolution
- Varied parameters: height, single/dual phase, partitioning




Energy Resolution



- Teflon reflectivity properties in GEANT4 class simulator
- PMTs with quantum efficiency of 35%
- 10% reflectivity grids
- 1-mm position resolution



Three rectangle configuration optimal at low energies (<300 keV)



Current Status

- Internals installed
- Outer can installed
- Xenon purified
- Gas in inner volume
- Detector cooling

Temp(K) of internals





- ✓ Cryogenic system installation and implementation
- ✓ Monte Carlo simulations demonstrating the energy resolution of a liquid xenon Compton imager
- Demonstration of energy resolution in PIXeY

Nicholas E. Destefano, Moshe Gai, Aaron G. Swindell, Alexander H. Young
University of Connecticut



Daniel N. McKinsey, Ethan P. Bernard, Sidney B. Cahn, Blair Edwards, Christopher J. Kachulis, Nicole A. Larsen, Alexey Lyashenko, James Nikkel, Christopher G. Wahl
Yale University



Cryogenics in PIXeY


- Impurities effect lifetime of drift electrons
- Decreases energy resolution
- < 1 ppb oxygen ideal

- Circulation system maintains Xenon purity
- LN₂ Cold Head provides cooling power
- Heat Exchanger reduces cooling power needed



✓ Has been installed, tested, and commissioned

Three Photon PET



Other 3 photon PET isotopes: ${}^{14}\text{O}$, ${}^{14}\text{C}$, ${}^{42\text{m}}\text{Sc}$, ${}^{48}\text{V}$, ${}^{50}\text{Mn}$, ${}^{52\text{m}}\text{Mn}$, ${}^{54}\text{Co}$, ${}^{55}\text{Co}$, ${}^{60}\text{Cu}$, ${}^{66}\text{Ga}$, ${}^{67}\text{Ge}$, ${}^{70}\text{As}$, and ${}^{73}\text{Se}$.

J.D. Kurfess *et al.*, IEEE, **2** (2001)1166.

Nicholas E. Destefano, UConn and Yale University. ANS&T Washington D.C. August 22, 2011