



Polymer-blend Organic Glass Scintillators for next generation neutron detectors at FRIB

Urmila Shirwadkar¹, Edgar V. van Loef¹, Tawan Jamdee¹, Lakshmi S. Pandian¹, Chuncheng Ji¹, Jarek Glodo¹, Patrick L. Feng², Annabelle Benin², Ryan Witzke², Remco G. T. Zegers³, Jorge Pereira³, and Cavan Maher³

¹*Radiation Monitoring Devices Inc., 44 Hunt Street, Watertown, MA 02472, USA*

²*Sandia National Laboratories, Livermore, CA 94550, USA*

³*Michigan State University, East Lansing, MI 48824, USA*

This work has been supported by the U.S. Department of Energy contract # DE-SC0021545. This support does not constitute an express or implied endorsement on the part of the Government. Distribution A. Approved for public release, distribution is unlimited.



Office of Science



Inspired by Light

rmd.dynasil.com

1

Radiation Monitoring Devices, Inc.

Mission

- Perform world class research
- Develop exceptional commercial products

Overview

- Founded 1974
- R&D and commercial products
- Acquired by Dynasil in 2008
- 80 employees

Research Expertise

- Semiconductors
- Scintillators
- Advanced Sensors
 - Nuclear
 - Optical
 - Magnetic
- Instruments & Systems
- Imaging
- Coatings



Organic Glass Scintillators for Nuclear Physics Experiments

Urmila Shirwadkar¹, Edgar V. van Loef¹, Tawan Jamdee¹, Lakshmi S. Pandian¹, Chuncheng Ji¹, Jarek Glodo¹, Patrick L. Feng², Annabelle Benin², Ryan Witzke², Remco G. T. Zegers³, Jorge Pereira³, and Cavan Maher³

¹*Radiation Monitoring Devices Inc., 44 Hunt Street, Watertown, MA 02472, USA*

²*Sandia National Laboratories, Livermore, CA 94550, USA*

³*Michigan State University, East Lansing, MI 48824, USA*

DOE Contract: DE-SC0021545 (Phase-I and II)

Administrative Contact: Michael Passo



Technical Contact: PM: Michelle Shinn

POC: Elizabeth Bartosz

POP: 4/4/2022 - 4/3/2024

Neutron detectors at FRIB

- Low Energy Neutron Detector Array (LENDA)* is used as a neutron time-of-flight (ToF) spectrometer developed for inverse kinematics (p,n) charge-exchange experiments in which low energy neutrons (~ 150 keV to 10 MeV) are detected.
- Versatile Array of Neutron Detectors at Low Energy (VANDLE) array is a hybrid neutron time-of-flight (TOF) detector consisting of plastic bars optimized for measurements of beta-delayed neutrons between 100 keV and 5 MeV
- LENDA Consists of 24 BC-408 plastic scintillator bars with dimensions 30 x 4.5 x 2.5 cm, each coupled to Hamamatsu H6410 PMT to each end of the bar
- These detector arrays are unable to differentiate between neutrons and gammas, complicating the background subtraction
- Separation between neutrons and photons is possible by time-of-flight, but it doesn't work for indirectly produced photons (e.g. produced by neutron-induced photons of surrounding materials)



LENDA+GRETINA@S800

picture: S. Noji

***Ref: G. Perdikakis et. al. LENDA: A low energy neutron detector array for experiments with radioactive beams in inverse kinematics, NIMA, Volume 686, 2012, Pages 117-124**

Ideal detector requirements

- Separation between neutrons and photons by pulse-shape discrimination (PSD)
- Need to cost-effectively cover large area
- Good timing resolution (< 1 ns)
- Position resolution (~ 5 cm)
- Low neutron-detection threshold (< 200 keV)



Issues with current detectors

- PSD slows down speed – affects timing
- Difficult to achieve with liquid scintillators that have PSD capability – toxicity, size, and mobility

Limited number of alternative options available and further development of capabilities is required

Organic glass scintillators as a solution

- Ability to provide PSD between gamma rays and neutrons
- Ability to scale up to large sizes
- High light yields of OGS allows for lower energy thresholds in PSD measurements and therefore a wider range of neutron energies to be detected
- Fast response
- Low cost






Scintillator	Emission (nm)	Light Yield (Ph/MeV)	Decay (ns)	PSD (FOM)
Stilbene	390	16,000	2.4-4.5	> 3
EJ-200	425	10,000	2.1	NO
EJ-232	370	8400	1.6	NO
EJ-232Q-1%	370	1700	0.7	NO
EJ-276	425	8,600	13	~ 2
Organic Glass Scintillators*	430	18,000	1.46	> 3

*Ref: Joseph S. Carlson, Patrick L. Feng, "Melt-cast organic glasses as high-efficiency fast neutron scintillators", Nucl. Instr. Meth. Phys. Res. A 832 (2016) 152-156

Tin-Loaded OGS

- Increased gamma-ray sensitivity, High Z_{eff} - spectroscopic capabilities
- Several reaction channels in one experiment, optimization of scientific output of complex and high-cost FRIB runs
- No compromise on the neutron detector solid angle to measure simultaneous gamma rays



Polymer-blended OGS

- Standard OGS are too brittle for high aspect ratio detectors that need to be self-supporting (30 x 4.5 x 2.5 cm³)
- Polymer blending with OGS can significantly enhance the mechanical properties of OGS*
- The core matrix material is an OGS-polymer mixture that comprises a small amount of polycarbonate or PVT blended into the OGS matrix.
- This polymer-blended OGS composition can be scaled up while maintaining high-performance similar to the standard OGS

****Ref: Nicholas R. Myllenbeck et. al. "Nano-segmented optical fibers containing molecular organic glass scintillator for fast neutron imaging", Proc. SPIE 11838, Hard X-Ray, Gamma-Ray, and Neutron Detector Physics XXIII, 118380R (1 September 2021); <https://doi.org/10.1117/12.2596532>***

Starting Materials, Synthesis and Purification

❖ Starting Materials

- ❑ Vinyltoluene/Polycarbonate - Alfa-Aesar/Fisher Scientific/Sigma-Aldrich
- ❑ OGS base material, bis(9,9'-dimethylfluoren-2-yl)diphenylsilane (P2) - Jubilant Biosys
- ❑ Organotin compound, trimethyl(4-vinylphenyl)stannane - Adesis, Inc. (New Castle, DE)
- ❑ bis-MSB - Sigma-Aldrich

❖ Synthesis

- ❑ Sandia prepared polycarbonate, Bis-MSB, and OGS compound
- ❑ RMD synthesized tin compound (TMS) as additives for incorporation into the OGS-Polymer blend


❖ Purification

- ❑ Done at RMD and Sandia

Polymer Organic Glass Scintillators


- ❑ Controlled polymerization of mixtures by dissolve OGS base material and a small amount of 1,4-bis[2-methylstyryl]-benzene (bis-MSB) into the styrene monomer or polycarbonate solution, and for tin version also add the organotin compound.
- ❑ Vary the concentrations of OGS and polymer to optimize polymer OGS composition
- ❑ Vary the concentrations of organotin compound to optimize light yield and PSD

Primary fluorophores




Diphenyl-bis(n,n-dimethylfluorene)silane

Secondary fluorophores




bis(methylstyryl)benzene

Polymer




Polystyrene

Polycarbonate doped OGS



AIB-6-101-B
1% PC, bis-MSB



AIB-6-116-2in
1% PC, BSBCz-Bu




BSO-PC3-2in
3% PC, bis-MSB


Sample	% PC	WLS	Rel LY EJ200	FOM
AIB-6-101-B	1	bis-MSB	1.40	2.31
AIB-6-116-2in	1	BSBCz-Bu	0.94	1.91
BSO-PC3-2in-PVA	3	bis-MSB	1.20	2.07

Large format Polymer OGS samples


3 wt% PC in OGS,
20 cm x 0.5 cm x 0.5 cm



1 wt% PC in OGS,
15 cm x 4.5 cm x 2.5 cm



3 wt% PS in OGS,
15 cm x 4.5 cm x 2.5 cm




Measurements at MSU

Objectives

- Evaluate minimum measurable energy-loss signal (threshold)
- Evaluate pulse-shape discrimination capabilities
- Evaluate timing resolution
- Measure neutron efficiency
- Compare results with LENDA bar



- DAQ used: NSCL/FRIB lab DAQ using two Pixie-16, 250 MHz digitizers
- ^{252}Cf fission source used for PSD and neutron efficiency measurements
- Various gamma sources used for detector characterization
- Well-developed ROOT-based analysis framework for LENDA detectors used*







Energy calibration and timing resolution

PSD and neutron efficiency measurements using two ELJEN liquid scintillators with the OGS samples

***Ref: S. Lipschutz et al., NIMA Volume 815, 1 (2016)**

Composition selection in Phase-I



Detector	Comments	Width of ^{241}Am 59.6 keV line (FWHM; keVee)	Low Threshold (in keVee)	Energy Limit	> 200 keVee	> 500 keVee	> threshold
RMD1	PC doped, 430 nm	17.59	~5 keV		1.22	1.51	.75
RMD2	PC doped, 530 nm	21.29	~20keV		.82	.88	.48
RMD3	Pure	14.22	~8 keV		1.08	1.31	.78
RMD4	Pure	13.69	~15 keV		1.23	1.27	N/A
RMD5	50% OGS in PVT	27.59	~28 keV		.62	.77	N/A
RMD6	Tin-loaded	19.59	~21 keV		1.14	1.185	N/A




Neutron efficiency measurement

- Efficiency determined relative to that of an ELJEN liquid scintillator (LSC) with known efficiency and normalized for detector volume
- Efficiency of RMD samples calculated by taking

$$Eff_{RMD}(E) = Eff_{LSCI} * \frac{N_{RMD}(E)}{N_{LSCI}(E)}$$



Measurements at MSU on 2" cylinder samples



AIB-6-101-B
1% PC, bis-MSB


AIB-6-116-2in
1% PC, BSBCz-Bu

BSO-PC3-2in
3% PC, bis-MSB


Detector	Sample number	Width of 59.6 keV line (FWHM; keV)	Low Threshold Limit (in keV)	22 Na First Compton Edge Channel #	PSD FOM > 200 keVee	PSD FOM: > 500 keVee
OGS-MS8	1	17.59	~5 keV	17924	1.22	1.51
OGS-3HF	2	21.29	~20keV	6837	.82	.88
SNL-170	3	14.22	~8 keV	16786	1.08	1.31
OGS-163.A2	4	13.69	~15 keV	17365	1.23	1.27
OGS-NRM-7-361.01	5	27.59	~28 keV	8872	.62	.77
OGS-TJ-008	6	19.59	~21 keV	9165	1.14	1.185
BSO-PC3	7	18.4	~30 keV	8772	0.98	1.02
AIB-6-101-B	8	16.38	~25 keV	13712	1.12	1.11
AIB-6-116	9	16.693	~30 keV	9256	1.22	1.18

Timing resolution




- Timing measurement performed relative to timing signal from a LENDA bar
- Resolution determined as a function of energy - contribution of LENDA timing measurement subtracted




Measurement of half bars



3% PS






1% PC





Measurement of half bars


3% PS



Dem: 0.492, FOM cut:[7, 50] Event#:172007



1% PS




Plans

- MSU to evaluate half bars
- Scale up of selected polymer OGS to 30 cm x 4.5 cm x 2.5 cm size
- Optimization of tin-loaded polymer OGS – use of vinyltoluene (aka methylstyrene)

Conclusions and Outlook

- LENDA is used for the detection of low-energy neutrons – options for upgrading the array with detectors that have pulse-shape discrimination for background reduction are being considered
- OGS composition has shown promise – scale up is possible due to polymer blending
- FRIB Charge-exchange group is characterizing half LENDA bar size PC-doped and PS-doped OGS samples are being evaluated at MSU
- Tin-loaded polymer OGS version is being developed

Acknowledgements



Work presented here is sponsored by DOE. This funding is greatly appreciated.

This work has been supported by the US Department of Energy contract # DE-SC0021545. This support does not constitute an express or implied endorsement on the part of the Government. Distribution A. Approved for public release, distribution is unlimited.